





**Dr. Maisch**  
Any Column, Any Size, Any Media



# CHIRAL COLUMNS

MADE BY DR. MAISCH



# CONTENT


- P 4      - REPROSIL CHIRAL POLYSACCHARIDE PHASES  
              - CHIRAL SELECTORS NAMES AND STRUCTURES
- P 5      - CROSS REFERENCE LIST  
              - IMMOBILIZED POLYSACCHARIDE PHASES
- P 6-11     - INDEPENDENT CHIRAL SFC STUDY BY SYAME KHATER  
              AND CAROLINE WEST, UNIVERSITY OF ORLEANS
- P 12-13    - COLUMN COMPARISON CHIRAL OM & CHIRACEL ODH
- P 14-15    - SELECTIVITY OF JM vs. OJ (USP-L80)
- P 16-17    - COMPARISON OF IMMOBILISED REPROSIL CHIRAL  
              POLYSACCHARIDE PHASES WITH COMPETITORS COATED  
              EQUIVALENTS
- P 17      - TYPICAL MOBILE PHASE CONDITIONS AND SOLVENT  
              COMPATIBILITY
- P 20-21    - REPROSIL CHIRAL NR AND NR-R
- P 22      - REPROSIL CHIRAL - AGP
- P 23-24    - PREPARATIVE COLUMNS
- P 25-30    - AVAILABLE COLUMNS AND PRODUCT NUMBERS

## CHIRAL COLUMNS MADE BY DR. MAISCH


From one of the biggest **High-Performance Liquid Chromatography (HPLC)** - Column Manufacturers in Europe.

# CHIRAL SELECTORS, NAMES & STRUCTURES OF COMMERCIALLY AVAILABLE CHIRAL STATIONARY PHASES (CSPs)

## „COATED“ REPROSIL CHIRAL POLYSACCHARIDE PHASES



## „IMMOBILISED“ REPROSIL CHIRAL POLYSACCHARIDE PHASES



## CROSS - REFERENCE IMMOBILIZED POLYSACCHARIDE PHASES


Dr Maisch phases	Daicel alternative	Phenomenex alternative	Phase description
ReproSil Chiral-MIA	Chiraldak IA	Lux i-Amylose-1	Amylose tris (3,5-dimethylphenylcarbamate)
ReproSil Chiral-MID	Chiraldak ID		Amylose tris (3-chlorophenylcarbamate)
ReproSil Chiral-MIF	Chiraldak IF		Amylose tris (3-chloro-4-methylphenylcarbamate)
ReproSil Chiral-MOF			Cellulose tris (4-chlorophenylcarbamate)
ReproSil Chiral-MIB	Chiraldak IB		Cellulose tris (3,5-dimethylphenylcarbamate)
ReproSil Chiral-MIC	Chiraldak IC	Lux i-Cellulose-5	Cellulose tris (3,5-dichlorophenylcarbamate)
ReproSil Chiral-MIX			Cellulose tris (4-chloro-3-methylphenylcarbamate)
ReproSil Chiral-MIZ			Cellulose (3-chloro-4-methylphenylcarbamate)
ReproSil Chiral-AM	Chiraldak AD	Lux Amylose-1	Amylose tris (3,5-dimethylphenylcarbamate)
ReproSil Chiral-AMS	Chiraldak AS		Amylose tris(S)-a-methylbenzylcarbamate)
ReproSil Chiral-YM	Chiraldak AY	Lux Amylose-2	Amylose tris (5-chloro-2-methylphenylcarbamate)
ReproSil Chiral-ZA	Chiraldak AZ		Amylose tris (3-chloro-4-methylphenylcarbamate)
ReproSil Chiral-BM	Chiralcel OB		Cellulose Tribenzoate
ReproSil Chiral-CM	Chiralcel OC		Cellulose tris (phenylcarbamate)
ReproSil Chiral-GM	Chiralcel OG		Cellulose tris (4-methylphenylcarbamate)
ReproSil Chiral-JM	Chiralcel OJ	Lux Cellulose-3	Cellulose tris (4-methylbenzoate)
ReproSil Chiral-OM	Chiralcel OD	Lux Cellulose-1	Cellulose tris (3,5-dimethylphenylcarbamate)
ReproSil Chiral-XM	Chiralcel OX	Lux Cellulose-4	Cellulose tris (4-chloro-3-methylphenylcarbamate)
ReproSil Chiral-ZM	Chiralcel OZ	Lux Cellulose-2	Cellulose tris (3-chloro-4-methylphenylcarbamate)

# MORE THAN 90% OF SEPARATIONS


REPROSIL CHIRAL COLUMNS ARE COMPLEMENTARY AND SUCCESSFUL IN ACHIEVING MORE THAN 90% OF SEPARATIONS OF THE ANALYTES.

ReproSil Chiral AM, JM and OM allow similar hit rates as competitor's CSPs.

## ReproSil CHIRAL COLUMNS



## DAICEL CHIRAL COLUMNS



- ReproSil Chiral AM
- ReproSil Chiral OM
- Non resolved enantiomers
- ReproSil Chiral JM

- Chiralkap AD
- Chiracel OD
- Non resolved enantiomers
- Chiracel OJ



# EXPERIMENTAL

## EVALUATION OF REPROSIL CHIRAL OM VS. OD

**Evaluation of a silica phase modified with cellulose tris-(3,5-dimethylphenyl- carbamate) „ReproSil Chiral-OM“ in supercritical fluid chromatography. Syame Khater and Caroline West, University of Orleans, CNRS UMR 7311, ICOA.**

All experiments were performed on a Jasco SFC system and an Acquity UPC<sup>2</sup> system. ReproSil Chiral-OM is based on silica coated with tris-(3,5-dimethylphenylcarbamate) of cellulose. Two hundred and thirty achiral compounds and one hundred and thirty chiral racemic compounds were screened on different polysaccharide-type chiral stationary phases in SFC in the following operating conditions: CO<sub>2</sub>/MeOH (90:10), flow rate 3 ml/min, oven temperature 25°C, outlet pressure 150 bars.


## NON-SPECIFIC INTERACTIONS AND RETENTION




Retention on cellulose tris-(3,5-dimethylphenylcarbamate) could be explained by non-specific interactions such as π-π interactions, hydrogen bonding and stereo-induced interactions.

# PERFORMANCE

Log k ReproSil OM





Log a CHIRALCEL OD



# APPLICATIONS

All experiments were performed on a Jasco SFC system and an Acquity UPC<sup>2</sup> system. ReproSil Chiral-OM is based on silica coated with tris-(3,5-dimethylphenylcarbamate) of cellulose. Two hundred and thirty achiral compounds and one hundred and thirty chiral racemic compounds were screened on different polysaccharide-type chiral stationary phases in SFC in the following operating conditions: CO<sub>2</sub>/MeOH (90:10), flow rate 3 ml/min, oven temperature 25°C, outlet pressure 150 bars.



The following chromatograms illustrate the complementary of the generic phases having cellulose tris-(3,5- dimethylphenyl-carbamate) as chiral selector in the course of method development: Focus on ReproSil Chiral-OM versus Chiralcel OD. The chromatograms illustrate the chiral compounds that are well resolved on ReproSil Chiral-OM (a) but have no separation on Chiralcel OD (b).



(2,3-epoxypropyl)-benzene

# IDENTICAL SELECTIVITY




## COMPARISON OF REPROSIL CHIRAL OM AND CHIRACEL OD






# COMPARISON OF USP-L40 CHIRAL COLUMNS

The analysis of fenamiphos, glutethimide and 5-methyl-5-phenylhydantoin on ReproSil Chiral-OM provide a better starting point for a method development than those on RegisCell, Lux Cellulose-1 or Cellucoat, respectively.

## NO SEPARATION





## SEPARATION



# COLUMN COMPARISON

## REPROSIL CHIRAL OM & CHIRACEL ODH



# REPROSIL CHIRAL OM & CHIRACEL ODH

## CONCLUSION

Reprosil Chiral OM is guaranteed replacement for Chiracel ODH HPLC columns.

## OBSERVATION

01.

The selectivity is equivalent either the columns used in basic, neutral or acidic conditions.

02.

Resolution between isomer is higher for ReproSil Chiral OM column when used in basic & neutral conditions where as in acidic condition Chiracel ODH is showing slightly higher resolution.

03.

The peak symmetry is better for ReproSil Chiral OM column in all three conditions (i.e. acidic,basic or neutral).



04.

With all the conditions ReproSil Chiral OM is showing higher no. of theoretical plate.



# SELECTIVITY OF JM vs. OJ (USP-L80)

The log  $\alpha$  - log  $\alpha$  plot compares the selectivities ability of JM and OJ.

The data points located on the first bisector (red dotted line) show similar separation profiles. The chromatograms of 1-(2-Naphthyl)-ethanol and Fonofos illustrate chiral compounds with identical separation factors on both columns.



The analysis of Hexobarbital on ReproSil Chiral-JM provides a better starting point for a method development than those on Chiral-AM or Chiral-OM.



## COMPARISON OF UNIQUE IMMOBILISED REPROSIL CHIRAL PHASES WITH COMPETITORS COATED EQUIVALENTS

### ReproSil Chiral – unique immobilised phases

#### SAMPLES

Sample: customer sample



Eluent: CO<sub>2</sub>/0 - 20% MeOH (0.1% DEA) in 2 min

#### ReproSil Chiral MIX

5 µm, 100 x 3.0 mm

#### Lux Cellulose-4

3 µm 100 x 4.6 mm



	RT	Areas	%Areas	Resolutions	USP Tailing
1	0,564	220320	39,72		1,20
2	0,726	334417	60,28	2,26	1,04

	RT	Areas	%Areas	Resolutions	USP Tailing
1	0,547	647055	44,58		
2	0,626	804384	55,42	0,97	



Sample: TSO  
Eluent: CO<sub>2</sub>/ 10% MeOH (0.1% DEA)

#### ReproSil Chiral MIX

5 µm, 100 x 3.0 mm

#### Lux Cellulose-4

3 µm 100 x 4.6 mm



	RT	Areas	%Areas	Resolutions	USP Tailing
1	0,588	1970619	49,64		1,13
2	0,789	1999432	50,36	4,97	1,11

	RT	Areas	%Areas	Resolutions	USP Tailing
1	0,583	1038369	49,91		1,05
2	0,851	1041945	50,09	8,89	1,04

## COMPARISON OF UNIQUE IMMOBILISED REPROSIL CHIRAL PHASES WITH COMPETITORS COATED EQUIVALENTS

### ReproSil Chiral – unique immobilised phases

#### SAMPLES

Sample: customer sample


Eluent: CO<sub>2</sub>/ 10-50% MeOH (0.1% DEA) in 2 min, hold until 5 min

#### ReproSil Chiral MIZ

3 µm, 100 x 3.0 mm


#### Lux Cellulose-2

3 µm 50 x 4.6 mm



	RT	Areas	%Areas	Resolutions	USP Tailing
1	1,585	1872263	36,75		1,06
2	1,699	3221889	36,25	1,77	1,09

	RT	Areas	%Areas	Resolutions	USP Tailing
1	1,678	1757077	36,89		1,05
2	1,801	3005764	63,11	2,00	1,06




[www.dr-maisch.com](http://www.dr-maisch.com)

## Comparison of unique immobilised ReproSil Chiral phases with competitor's coated equivalents in NP and PO mode

### ReproSil Chiral MIX vs. Lux Cellulose-4 under NP/PO conditions

Customer sample  
No additives used  
Black: Lux Cellulose 4  
**Colour: ReproSil Chiral MIX**



## REPROSIL CHIRAL – TYPICAL MOBILE PHASE CONDITIONS

### Normal Phase

- Mixtures of hexane or heptane with alcohols (EtOH, IPA) = 80:20 (vary % alcohol to adjust retention time and selectivity)
- Add 0.1 – 0.5% DEA or TEA for basic analytes and 0.1 – 0.5 % TFA or AcOH for acidic analytes

### Polar Organic Phase

- Mixtures of ACN / IPA (95/5) or MeOH / IPA (90/10) or neat ACN
- Add 0.1 – 0.5% DEA or TEA for basic analytes and 0.1 – 0.5 % TFA or AcOH for acidic analytes

### Reversed Phase

- ACN or MeOH or EtOH / water mixtures
- Water content must be < 85%
- Add 0.5 – 1 N perchlorate or 0.1% TFA for basic compounds and HClO<sub>4</sub>/NaClO<sub>4</sub> buffer for acidic compounds together with ACN
- Use dedicated column for reversed phase conditions because solvent switch is tedious

### Solvent compatibility of coated phases

As the polysaccharide layer is not chemically bonded to the silica support but physically adsorbed, strong solvents have the ability to strip the polymer off the silica support and have therefore to be avoided even in trace amounts. Such strong solvents include:

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>Ethers incl. THF</li> <li>Acetone</li> <li>Chlorinated solvents</li> <li>Ethyl acetate</li> <li>DMSO</li> </ul> | <ul style="list-style-type: none"> <li>DMF</li> <li>Toluene</li> <li>Ketones</li> <li>Dimethylacetamid</li> <li>IPA &gt; 50%</li> </ul> |
|--|---|

**With immobilised phases an extended range of organic solvents can be used as injection solvents and eluents.**

Recommended temperature range: 0 – 40°C

Recommended max. pressure: 150 bar


# ReproSil Chiral NR – unique chemically bonded phase

- Immobilised brush-type phase
- Dinitro compound, acceptor and electron donor phase
- Particularly useful for aromatic compounds with O or N near chiral centre
- NR, RP, SFC mode
- Both antipodes of chiral selector available - elution order reversible

## Why ReproSil Chiral NR?

- switch from **NP** to **RP** within minutes
- In all common HPLC eluents chemically stable
- amenable to various separation modes and solvents
- scalable
- complementary selectivity
- elution order reversible




## CHIRAL PREPARATIVE SEPARATIONS WITH HPLC BY HOFFMANN LA ROCHE, BASEL IN 2014



## Comparison of ReproSil Chiral NR vs. Whelk-O-1 in SFC mode

### ReproSil Chiral NR vs. Whelk-O1

COMPARISON



Co-Solvent:  
MeOH(0.1%DEA)

Gradient (B%):  
10% to 50% in 2.0 min.  
hold 1.0 min at 50%

Temperature: 35°C

Flow (ml/min): 4

Back Pressure(psi):  
1500.00

Detector: 220nm

# CHIRAL SEPARATION OF IBUPROFEN, WARFARIN AND VERAPAMIL




# LONGLIFE® Preparative Chromatography Column Hardware

ReproSil Chiral-AGP, 5 µm

ELUENTS

Technology patented by  
Dr. Maisch HPLC

LONGLIFE




Ibuprofen  
100 mM  
Na-phosphate buffer, pH7

Warfarin  
10% -Propanol in 10 mM  
Na-phosphate buffer, pH7

Verapamil  
12% ACN in 10 mM  
Na-phosphate buffer, pH7

THE ONLY PREPACKED COLUMN ON THE MARKET  
WITH INTEGRATED DYNAMIC AXIAL  
COMPRESSION MECHANISM

Available column IDs: 25 mm, 30 mm, 40 mm, 50 mm, 70 mm







[www.dr-maisch.com](http://www.dr-maisch.com)

# Preparative columns Quality Assurance Chromatogram

## CHIRAL COLUMNS

### ReproSil Chiral-NR, 8 µm upscaling



#### TEST CONDITIONS

Mobile Phase:  
Heptane/IPA 85/15  
Temperature: Ambient  
Pressure: 15 bar  
Detector: UV @ 229 nm  
Sensitivity: 0,1 mV

Mobile Phase:  
Heptane/IPA 85/15  
Temperature: Ambient  
Pressure: 9 bar  
Detector: UV @ 254 nm  
Sensitivity: 0,0 mV

Mobile Phase:  
Heptane/IPA 85/15  
Temperature: Ambient  
Pressure: 11 bar  
Detector: UV @ 254 nm  
Sensitivity: -0,8 mV

Mobile Phase:  
Heptane/IPA 85/15  
Temperature: Ambient  
Pressure: 12 bar  
Detector: UV @ 254 nm  
Sensitivity: -0,3 mV

Exemplary columns are mentioned below.  
Columns with ID of: 4.0 mm, 3.0 mm and 2.0 mm are also available.

#### Immobilized Amylose-Phases

<b>ReproSil Chiral-MIA</b>	<b>Amylose tris-(3,5-dimethylphenyl) carbamate</b>	(PN: r65.mia.)
Chiral-MIA, 5µm (or ReproSil Chiral-MIA-R)	(Alternative to Daicel IA) 250 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm

<b>ReproSil Chiral-MIF</b>	<b>Amylose tris-(3-Chloro-4-Methylphenyl) carbamate</b>	(PN: r65.mif.)
Chiral-MIF, 5 µm (or ReproSil Chiral-MIF-R)	(Alternative to Daicel IF) 250 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm

<b>ReproSil Chiral-MID</b>	<b>Amylose tris-(3-Chlorophenyl) carbamate</b>	(PN: r65.mid.)
Chiral-MID, 5 µm (or ReproSil Chiral-MID-R)	(Alternative to Daicel ID) 250 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm

#### Immobilized Cellulose-Phases

<b>ReproSil Chiral-MIB</b>	<b>Cellulose tris-(3,5-dimethylphenyl) carbamate</b>	(PN: r65.mib.)
ReproSil Chiral-MIB, 5 µm (or ReproSil Chiral-MIB-R)	(Alternative to Daicel IB) 250 x 4.6 mm	150 x 4.6 mm

<b>ReproSil Chiral-MIC</b>	<b>Cellulose tris-(3,5-Dichlorophenyl) carbamate</b>	(PN: r65.mic.)
ReproSil Chiral-MIC, 5 µm (or ReproSil Chiral-MIC-R)	(Alternative to Daicel IC) 250 x 4.6 mm	150 x 4.6 mm
ReproSil Chiral-MIC, 3 µm (or ReproSil Chiral-MIC-R)	250 x 10 mm (Alternative to Daicel IC-3) 250 x 4.6 mm	250 x 20 mm 150 x 4.6 mm

<b>ReproSil Chiral-MIX</b>	<b>Cellulose tris-(4-Chlor-3-Methylphenyl) carbamate</b>	(PN: r65.mix.)
ReproSil Chiral-MIC, 5 µm (or ReproSil Chiral-MIX-R)	250 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm

<b>ReproSil Chiral-MIZ</b>	<b>Cellulose tris-(3-Chlor-4-methylphenyl) carbamate</b>	(PN: r65.miz.)
ReproSil Chiral-MIZ, 5 µm (or ReproSil Chiral-MIZ-R)	250 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm
ReproSil Chiral-MIZ, 3 µm (or ReproSil Chiral-MIZ-R)	250 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm

# CHIRAL COLUMNS

<b>ReproSil Chiral-MOF</b>	<b>Cellulose tris-(4-Chlorphenyl) carbamate</b>		
ReproSil Chiral-MOF, 5 µm (or ReproSil Chiral-MOF-R)	250 x 4.6 mm	150 x 4.6 mm	(PN: r65.mof.) 100 x 4.6 mm
<b>Coated Cellulose Phases</b>			
<b>ReproSil Chiral-OM</b>	<b>USP-L40, Cellulose tris-(3,5-dimethylphenyl-carbamate) mod. Silica</b>		
ReproSil Chiral-OM, 3 µm / ReproSil Chiral-OM-R, 3 µm	Alternatives to Daicel OD-3 / OD-3R		
250 x 4.6 mm 100 x 4.6 mm	150 x 4.6 mm 50 x 4.6 mm	125 x 4.6 mm 33 x 4.6 mm	
ReproSil Chiral-OM, 5 µm / ReproSil Chiral-OM-R, 5 µm	Alternative to Daicel OD-H / OD-RH		
250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 250 x 20 mm	100 x 4.6 mm	
ReproSil Chiral-OM, 10 µm / ReproSil Chiral-OM-R, 10 µm	Alternatives to Daicel OD / OD-R		
250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 250 x 20 mm	100 x 4.6 mm	
ReproSil Chiral-OM, 20 µm / ReproSil Chiral-OMR, 20 µm	Alternative to Daicel OD / OD-R		
250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 250 x 20 mm	100 x 4.6 mm	
<b>ReproSil Chiral-CM</b>	<b>USP-L70 (Tris-(Phenylcarbamate)-Cellulose mod. Silica)</b>		
ReproSil Chiral-CM, 3 µm	Alternative to Daicel OC-3 250 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm	(PN: r63.cm)
ReproSil Chiral-CM, 5 µm	Alternative to Daicel OC-H 250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 100 x 4.6 mm	(PN: r65.cm)
<b>ReproSil Chiral-JM</b>	<b>USP-L80 (Tris-(4-Methylbenzoyl)-Cellulose mod. Silica)</b>		
ReproSil Chiral-JM, 3 µm (or ReproSil Chiral-JM-R)	Alternative to Daicel OJ-3 / OJ-3R 250 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm	(PN: r63.jm)
ReproSil Chiral-JM, 5 µm (or ReproSil Chiral-JM-R)	Alternative to Daicel OJ-H / OJ-HR 250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 100 x 4.6 mm	(PN: r65.jm)
ReproSil Chiral-JM, 10 µm (or ReproSil Chiral-JM-R)	Alternative to Daicel OJ / OJR 250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 100 x 4.6 mm	(PN: r60.jm)

# CHIRAL COLUMNS

<b>ReproSil Chiral- ZM</b>	<b>Cellulose tris-3-Chloro-4-Methylphenylcarbamate mod. Silica</b>		
ReproSil Chiral-ZM, 3 µm (or ReproSil Chiral-ZM-R)	Alternative to Daicel OZ-3 / OZ-3H 250 x 4.6 mm 100 x 4.6 mm	150 x 4.6 mm 125 x 4.6 mm 50 x 4.6 mm	(PN: r63.zm) 100 x 4.6 mm 250 x 20 mm
<b>Coated Amylose Phases</b>			
<b>ReproSil Chiral-BM</b>	<b>Tris-(Benzoyl)-Cellulose mod. Silica</b>		
ReproSil Chiral-BM, 5 µm / ReproSil Chiral-BM-R	Alternatives to Daicel OB-H / OB-RH 250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 100 x 4.6 mm	(PN: r65.zm)
ReproSil Chiral-BM, 10 µm / ReproSil Chiral-BM-R	Alternatives to Daicel OB / OB-R 250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 100 x 4.6 mm	150 x 4.6 mm 100 x 4.6 mm
<b>ReproSil Chiral-AM</b>	<b>USP-L51, Amylose tris-3,5-dimethylphenylcarbamate mod. Silica</b>		
ReproSil Chiral-AM, 3 µm / ReproSil Chiral-AM-R, 3 µm	Alternative to Daicel AD-3 / AD-3R 250 x 4.6 mm 100 x 4.6 mm	150 x 4.6 mm 50 x 4.6 mm	125 x 4.6 mm 33 x 4.6 mm
ReproSil Chiral-AM, 5 µm / ReproSil Chiral-AM-R, 5 µm	Alternatives to Daicel AD-H / AD-RH 250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 250 x 20 mm	100 x 4.6 mm
ReproSil Chiral-AM, 10 µm / ReproSil Chiral-AM-R, 10 µm	Alternatives to Daicel AD / AD-R 250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 250 x 20 mm	100 x 4.6 mm
ReproSil Chiral-AM, 20 µm / ReproSil Chiral-AM-R, 20 µm	Alternatives to Daicel AD / AD-R 250 x 4.6 mm 250 x 10 mm	150 x 4.6 mm 250 x 20 mm	100 x 4.6 mm

# CHIRAL COLUMNS

## **ReproSil Chiral-AMS**

### **Amylose tris-(S)- $\alpha$ -Methylbenzyl-Carbamate**

ReproSil Chiral-AMS, 3 $\mu$ m / ReproSil Chiral-AMS-R, 3 $\mu$ m	Alternatives to Daicel AS-3 / AS-3R
250 x 4.6 mm	150 x 4.6 mm
100 x 4.6 mm	50 x 4.6 mm
ReproSil Chiral-AMS, 5 $\mu$ m / ReproSil Chiral-AMS-R, 5 $\mu$ m	Alternatives to Daicel AS-H / AS-RH
250 x 4.6 mm	150 x 4.6 mm
250 x 10 mm	100 x 4.6 mm
ReproSil Chiral-AMS, 10 $\mu$ m / ReproSil Chiral-AMS-R, 10 $\mu$ m	Alternatives to Daicel AS / AS-R
250 x 4.6 mm	150 x 4.6 mm
250 x 10 mm	100 x 4.6 mm

## **ReproSil Chiral-ZA**


### **Amylose tris-3-Chloro-4-Methylphenylcarbamate mod. Silica**

Alternative to Daicel AZ-3 / AZ-3R	(PN: r63.za.)
250 x 4.6 mm	150 x 4.6 mm
100 x 4.6 mm	125 x 4.6 mm
	50 x 4.6 mm
	33 x 4.6 mm
Alternative to Daicel AZ-H / AZ-HR	(PN: r65.za.)
250 x 4.6 mm	150 x 4.6 mm
250 x 10 mm	100 x 4.6 mm
	250 x 20 mm

## **ReproSil Chiral-YM**

### **Amylose tris-5-Chloro-2-Methylphenylcarbamate mod. Silica**

Alternatives to Daicel AY-3 / AY-3R	(PN: r63.ym.)
250 x 4.6 mm	150 x 4.6 mm
100 x 4.6 mm	125 x 4.6 mm
	50 x 4.6 mm
	33 x 4.6 mm



# CHIRAL COLUMNS

## **ReproSil Chiral-NR**

Immobilized brush-type phases, Dinitro-compounds.  
 $\pi$ -electron acceptor /  $\pi$ -electron donor phase. Particularly for aromatic compounds with O or N near chiral-centre. Chiral separations in NP and RP-Modus

5 $\mu$ m	(PN: r15.nr)	250 x 4.6 mm
		150 x 4.6 mm
		100 x 4.6 mm
8 $\mu$ m	(PN: r18.nr)	250 x 4.6 mm
12 $\mu$ m	(PN: r112.nr)	150 x 4.6 mm
15 $\mu$ m	(PN: r115.nr)	100 x 4.6 mm
		250 x 10 mm
		250 x 20 mm

## **ReproSil Chiral-NR-R**

Reversed Elution order compared to ReproSil Chiral-NR / Antipode of Chiral-NR

(PN: r18.nr)	250 x 4.6 mm
(PN: r112.nr)	150 x 4.6 mm
	100 x 4.6 mm
	250 x 10 mm
	250 x 20 mm

## **Other Chiral Phases**

### **ReproSil Chiral-PS, 8 $\mu$ m**

Chiral separations in NP and RP-Modus  
 (for aromatic compounds with S or P, for example: Sulfoxides, Phosphine Oxides, Phosphonates, Thiophosphine Oxides, Phosphin Selenides, Phosphine Boranes)

250 x 4.0 mm	(PN: r18.ps.)
--------------	---------------

### **ReproSil Chiral-OH, 8 $\mu$ m**

Chiral separations in NP-Modus  
 (Aromatic Alcohols with OH near chiral centre, Aryl Carbinols)

250 x 4.0 mm	(PN: r18.oh.)
--------------	---------------

### **ReproSil Chiral-AA, 8 $\mu$ m**

(for all amino acids, L-form elutes first)

250 x 4.0 mm	(PN: r18.aa.)
--------------	---------------

### **ReproSil Chiral-TAG, 8 $\mu$ m**

Chiral separations in NP and RP-Modus

(Teicoplanin Aglycon: for Aminoalkohols, N-blocked amino acids, a-Hydroxy Acids, Oxazolidinons, Hydantoins, Imides, Amino Acids USP-L63)

250 x 4.0 mm	(PN: r18.tag.)
--------------	----------------

### **ReproSil Chiral-Beta-CD, 5 $\mu$ m**

USP-L45 (Dansyl-Aminoacids, Barbiturates, Propranolol Sulfonamide, Prostaglandines)

# CHIRAL COLUMNS

**ReproSil Chiral-Gamma-CD, 5 µm** 250 x 4.0 mm (PN: r15.gcd.s2504)

**ReproSil Chiral-D-PhenylGlycin, 5 µm** USP-L36, N-(3,5-Dinitrobenzoyl)-D-Phenylglycin.

(Herbicides + Pharmaca (Alcohols, Carbon. acids, Esters, Sulfoxides) Fenoprop-Methyl, Mecoprop-Methyl, Supidimid)

250 x 4.0 mm (PN: r15.DPG.s2504)  
250 x 10 mm  
250 x 20 mm

**ReproSil Chiral-L-PhenylGlycin, 5 µm** USP-L36 N-(3,5-Dinitrobenzoyl)-L-Phenylglycin.  
Antipode to D-PhenylGlycin (Herbicides + Pharmaca (Alcohols, Carbon. Acids, Esters, Sulfoxides)  
Fenoprop-Methyl, Mecoprop-Methyl, Supidimid)

250 x 4.0 mm (PN: r15.LPG.s2504)

**ReproSil Chiral-L-Leucin, 5 µm** N-(3,5-Dinitrobenzoyl)-L-Leucin  
250 x 4.0 mm (PN: r15\_LL.s2546)


**ReproSil Chiral-L-Prolin, 5 µm** Davankov-Ligand exchange, Aminoacids, Hydantoin, Succinimide, Gluthetimide, Barbiturate, Sulfoxide  
250 x 4.0 mm (PN: r15.pr.s2504)

**ReproSil L-Hydroxy-Prolin, 5 µm** Alternative to Nucleosil Chiral-1, USP-L32  
Ligand exchange, Eluent: 2-10 mM Coppersulfate, 20-60 C° (For DL-Atrolactinacid, DL-Mandelic acids, DL-Lactic acid, DL-Asparagin, DL-Serin, DL-Phenylalanin, DL-Threonin, DL-Prolin, DL-Histidin, DL-Valin, DL-Tyrosin, DL-Tryptophan),  
250 x 4.0 mm (PN: r15.hp.s2504)

**ReproSil Chiral-AGP, 5 µm** 300 Å Silicagel with chiral AGP-Protein, USP-L41

100 x 4.0 mm	(PN: r35.agp.s1004)	150 x 2.0 mm	(PN: r35.agp.s1502)
150 x 3.0 mm	(PN: r35.agp.s1503)	100 x 2.0 mm	(PN: r35.agp.s1002)
100 x 3.0 mm	(PN: r35.agp.s1003)	50 x 2.0 mm	(PN: r35.agp.s0502)
50 x 3.0 mm	(PN: r35.agp.s0503)		
Guards: 5 x 3.0 mm: 2 pieces	(PN: r35.agp.v0003)	Guards: 5 x 2.0 mm: 2 pieces	(PN: r35.agp.v0002)
Guard-holder, direct	(PN: 81.00)	Guard-holder, direct	(PN: 91.00)

# CHIRAL COLUMNS



# **Dr. Maisch**

Any Column, Any Size, Any Media

Distributor:

Dr. Maisch HPLC GmbH  
Beim Brückle 14  
D-72119 Ammerbuch  
T: +49 7073 50357  
F: +49 7073 4216  
[www.Dr-Maisch.com](http://www.Dr-Maisch.com)  
[www.MODcol.com](http://www.MODcol.com)  
E-Mail: [info@dr-maisch.com](mailto:info@dr-maisch.com)



PDF brochure  
for download

