

Dr. Maisch

Any Column, Any Size, Any Media

ReproSil XR
eXtRa High Purity Phases



ReproSil™ XR
HPLC Columns

- Highly Efficient and Perfect Reproducibility
- General Use Selectivity
- Exceptional Value

ReproSil™ XR	Specifications
Phase	Si, Amino, Phenyl, C8 and C18 (MS)
Encapping	yes
Typical Efficiency	>90,000 plates/m
Surface Area	250 m ² /g , 100 m ² /g
Pore Size	120Å, 300Å
Particle Size	1.5 ,3 ,3.5 ,5 ,7 and 10µm

Choose Wisely - ReproSil XR

Dr.Maisch HPLC GmbH sets a new standard in HPLC columns. We make everything - from the silica particle to the finished product. Our streamlined manufacturing process allows delivery of the highest performance at an exceptional value!

ReproSil XR Columns			
Part No.	i.d. x Length	Particle Size	Price
rx15.9e.s1546	4.6 x 150mm	5µm	240€
rx15.9e.s2546	4.6 x 250mm	5µm	270€
rx15.9e.s1502	2.0 x 150mm	5µm	210€
rx15.9e.s2502	2.0 x 250mm	5µm	260€


XtRa High Quality HPLC Phases

Extra High purity silica and uniform bonded phase coverage translates to symmetrical peaks for acids/bases, and predictable reversed phase selectivity. Whether routine analysis or new method developments, use ReproSil XR columns to get premium performance without a premium price.

Reproducible Methods Start with Reproducible Columns


Our tightly controlled silica synthesis and bonding keep capacity factor and selectivity variations to a minimum. The advanced packing methods deliver consistently high column-to-column efficiency.

Highly Reproducible Selectivity



Consistent batch-to-batch capacity for acidic, basic & neutral compounds

Standards Mix




1. Nicotine
2. Quinine
3. Diphenhydramine
4. Nortriptyline
5. 2,4-Dichlorophenoxy acetic acid
6. 2,4-Dichlorophenoxy propionic acid
7. 2,4-Dichlorophenoxy butyric acid

Column: ReproSil XR 120 C18
5µm, 150 x 4.6mm
MobilePhase: 50% MeCN / 50% 50mM
KH₂PO₄ pH 3 at 35°C

Baseline resolve challenging separations

Acids, Bases and Neutrals




1. Uracil
2. Amitriptyline
3. Phenylvaleric Acid
4. Toluene

Column: ReproSil XR 120 C18,
5µm, 150 x 4.6mm
MobilePhase: 50% MeCN/50% 50mM
KH₂PO₄ pH 3 at 35°C

Get symmetrical peaks shapes for strong acids and bases


Expect Competitive Performance

Efficient manufacturing is driving force for competitive pricing, but certainly not at the risk of non-competitive performance. ReproSil XR columns show equal or better efficiency and asymmetries for challenging base and acid components in comparison to industry leading columns.



ReproSil XR columns compares favorable to leading competitor columns

Amoxicillin (U.S.P.)




— ReproSil XR C18 Column
 — Competitor P C18 Column
 — Competitor G C18 Column
 — Competitor H C18 Column
 — Competitor K C18 Column

Column: ReproSil XR C18 Column, 4.6 x 250mm, 5µm
 (Part No. **rx15.9e.s2546**)
MobilePhase: Isocratic, 0.05 N monobasic potassium phosphate pH 5:
 Acetonitrile 96:4
FlowRate: 1.5ml/min
Detector: UV, 230nm
ColumnTemp: 30°C
InjectionVol: 10µl


	USP Spec	ReproSil XR C18 Col	Comp P C18 Col	Comp G C18 Col	Comp H C18 Col	Comp K C18 Col
Ret. Time (min)	N/A	3.80	4.05	5.09	3.34	3.86
Theoretical plates	≥1700	8650	7240	6230	7060	7680
Tailing Factor	≤2.5	1.0	1.2	1.1	1.0	1.0
k'	1.1≤k'≤2.8	1.2	1.6	2.3	1.0	1.6

Water Soluble Vitamins



Column: ReproSil XR 120 C18 Column, 4.6 x 150mm
 (Part No. **rx15.9e.s1546**)
MobilePhase: A: 0.1M NH₄H₂PO₄ Buffer, pH4.5
 B: Acetonitrile
Gradient: (Time, %B): (0,5), (4,5), (10,40), (15,40)
FlowRate: 1.0mL/min
Detector: UV at 254nm
InjectionVol: 2µl


Phenols



Column: ReproSil XR 120 C18 Column, 4.6 x 150mm
 (Part No. **rx15.9e.s1546**)
MobilePhase: A: 1% Glacial Acetic Acid in Water
 B: 1% Glacial Acetic Acid in Methanol
Gradient: (Time, %B): (0,45), (6,45) (8,60), (15,80)
FlowRate: 1.0mL/min
Detector: UV at 280nm
InjectionVol: 10µl


Universally usable!

Organic Acids




Column: ReproSil XR 120 C18 Column, 120Å, 5µm, 4.6 x 150mm
 (Part No. **rx15.9e.s1546**)
MobilePhase: 0.025M KH₂PO₄(pH2.5) : Methanol, 97:3
FlowRate: 0.7ml/min
Detector: UV at 220nm
InjectionVol: 1µl

Carbamate & Urea




Column: ReproSil XR 120 C18 Column, 120Å, 5µm, 4.6 x 150mm
 (Part No. **rx15.9e.s1546**)
MobilePhase: A: 0.025M KH₂PO₄, pH3.2 B: Acetonitrile
Gradient: (Time, %B): (0,30), (20,60)
FlowRate: 1.0ml/min
Detector: UV at 240nm
InjectionVol: 10µl

Paracetamol Effervescent Tablets on Chinese Pharmacopeia Assay




Column: ReproSil XR 120 C18, 5µm, 250x4.6mm
Mobile Phase: Phosphate Buffer:Methanol(80:20)
Flow Rate: 1.0mL/min
Detector: UV at 254nm

Folium Ginkgo




Column: ReproSil XR 120 C18, 5µm, 250x4.6mm
Mobile Phase: Methanol:Tetrahydrofuran:Water(25:10:65)
Flow Rate: 1.0mL/min
Detector: ELSD

Simvastatin Tablets (U.S.P)



Column: ReproSil XR 120 C18 Column, 4.6 x 250mm
 (Part No. **rx15.9e.s2546**)
MobilePhase: 0.39% monobasic sodium phosphate pH4.5: Acetonitrile 35:65
FlowRate: 1.5ml/min
Detector: UV, 238nm
ColumnTemp: 45°C
InjectionVol: 10µl

Ibuprofen (U.S.P)



Column: ReproSil XR 120 C18 Column, 4.6 x 150mm
 (Part No. **rx15.9e.s1546**)
MobilePhase: H2O (adjusted with phosphoric acid to a pH of 2.5) :
 Acetonitrile 6.5:3.5
FlowRate: 2ml/min
Detector: UV, 214nm
InjectionVol: 5µl

From UHPLC

to PrepLC

fully scaleable



The information presented herein is derived from our testing and experience. It is offered, free of charge, for your consideration, investigation and verification. Since operating conditions vary significantly, and since they are not under our control, we disclaim any and all warranties on the results which might be obtained from the use of our products. You should make no assumption that all safety or environmental protection measures are indicated or that other measures may not be required. This product(s) may be covered by patents or patents pending. Dr.Maisch HPLC GmbH reserves the right to change prices and/or specifications without prior notification. Pictures licensed by Shutterstock.com. Printed in the Germany.

Dr. Maisch
 Any Column, Any Size, Any Media

Distributor:

Dr. Maisch HPLC GmbH
 Beim Brückle 14
 D-72119 Ammerbuch
 Fon: 0049 7073 50357
 Fax: 0049 7073 4216
www.Dr-Maisch.com
www.ModCol.com
 E-Mail: info@Dr-Maisch.com



PDF brochure
 for download