

Column comparison sheet

Manufacturer	TYPE	Column	Direct replacement	1st alternative	2nd alternative
PHENOMENEX	Core-Shell	AERIS WIDEPORE C4	n/a	n/a	Reprosil-XR 300 C4
PHENOMENEX	Core-Shell	AERIS WIDEPORE XB-C8	n/a	n/a	Reprosil-XR 300 C4
PHENOMENEX	Core-Shell	AERIS WIDEPORE XB-C18	n/a	n/a	Reprosil-XR 300 C18
PHENOMENEX	Core-Shell	AERIS PEPTIDE XB-C18	n/a	n/a	Reproshell ODS-1
PHENOMENEX	Fully-porous silica	AQUA 125 C18	n/a	Reprospher Aqua C18	n/a
PHENOMENEX	Fully-porous silica	AQUA 200 C18	n/a	Reprospher Aqua C18	n/a
PHENOMENEX	Fully-porous silica	BIOSEP 2000	ReproSil 200 SEC-2	Reprosil 200 SEC	n/a
PHENOMENEX	Fully-porous silica	BIOSEP 3000	Reprosil 300 SEC	n/a	n/a
PHENOMENEX	Fully-porous silica	BIOSEP 4000	Reprosil 4000 SEC	n/a	n/a
PHENOMENEX	Fully-porous silica	BIOZEN SEC-2	n/a	ReproSil 200 SEC-2	Reprosil 200 SEC
PHENOMENEX	Fully-porous silica	BIOZEN SEC-3	n/a	Reprosil 300 SEC	n/a
PHENOMENEX	Fully-porous silica	BIOZEN Glycan	n/a	n/a	n/a
PHENOMENEX	Fully-porous silica	BIOZEN Intact C4	n/a	n/a	Reprosil-XR 300 C4
PHENOMENEX	Fully-porous silica	BIOZEN Intact XB-C8	n/a	n/a	Reprosil-XR 300 C4
PHENOMENEX	Fully-porous silica	BIOZEN Peptide PS-C18	n/a	n/a	n/a
PHENOMENEX	Core-Shell	BIOZEN Peptide XB-C18	n/a	n/a	n/a
PHENOMENEX	Non-porous polymer	BIOZEN WCX	n/a	n/a	n/a
PHENOMENEX	Fully-porous silica	BONDCLONE C18	Reprobond C18	n/a	n/a
PHENOMENEX	Fully-porous silica	BONDCLONE Silica	Reprobond Si	n/a	n/a
PHENOMENEX	Fully-porous silica	CHIREX 3014	n/a	n/a	n/a
PHENOMENEX	Core-Shell	CLARITY OLIGO-MS	n/a	Reproshell ODS-1	n/a
PHENOMENEX	Ethylene-bridged silica	CLARITY Oligo-RP	Reprosil Phoenix C18	n/a	n/a
PHENOMENEX	NPS Polymer	CLARITY Oligo-SAX	n/a	n/a	n/a
PHENOMENEX	Core-Shell	CLARITY Oligo-XT	n/a	n/a	Reproshell ODS-1
PHENOMENEX	CHIRAL-Pirkle brush	CHIREX 3020	n/a	n/a	n/a
PHENOMENEX	Dankov phase	CHIREX 3126	n/a	n/a	Reprosil L-Hydroxy-Prolin,
PHENOMENEX	CHIRAL-Pirkle brush	CHIREX 3001	n/a	n/a	n/a
PHENOMENEX	CHIRAL-Pirkle brush	CHIREX 3019	n/a	n/a	n/a
PHENOMENEX	CHIRAL-Pirkle brush	CHIREX 3022	n/a	n/a	n/a
PHENOMENEX	Fully-porous silica	COLUMBUS C8	n/a	Reprospher 100 C8-TDE	n/a
PHENOMENEX	Fully-porous silica	COLUMBUS C18	n/a	Reprospher 100 C18-TDE	n/a

Dr. Maisch LC column replaces Phenomenex column 1:1. Small adaptations of the existing method may be necessary.

Dr. Maisch LC column could replace LC column. Adaptation of the method probably necessary.

Technical specifications of the Dr. Maisch LC column are the closest to the Phenomenex column. Small risk that column does not represent a good replacement.

PHENOMENEX	Fully-porous silica	CUROSIL	n/a	ReproSil Fluosil 100 PFP	n/a
PHENOMENEX	Fully-porous silica	ENVIROSEP PP PAH	n/a	Ultrasep ES PAH Reprosil PAH	
PHENOMENEX	Ethylene-bridged silica	GEMINI C18	Reprosil Phoenix C18	Reprosil Gold 120 C18	n/a
PHENOMENEX	Ethylene-bridged silica	GEMINI C6-Phenyl	n/a	n/a	Reprospher 100 Phenyl-hexyl-e
PHENOMENEX	Ethylene-bridged silica	GEMINI-NX C18	Reprosil Phoenix C18	Reprosil Gold 120 C18	n/a
PHENOMENEX	Fully-porous silica	HYPERCLONE BDS C18	Equisil BDS C18	n/a	n/a
PHENOMENEX	Fully-porous silica	HYPERCLONE BDS C8	Equisil BDS C8	n/a	n/a
PHENOMENEX	Fully-porous silica	HYPERCLONE CN (CPS)	Equisil CPS	n/a	n/a
PHENOMENEX	Fully-porous silica	HYPERCLONE MOS (C8)	Equisil MOS-2	n/a	n/a
PHENOMENEX	Fully-porous silica	HYPERCLONE ODS (C18)	Equisil ODS	n/a	n/a
PHENOMENEX	Fully-porous silica	INERTCLONE C8	n/a	Reprosil Pur 120 C8	
PHENOMENEX	Fully-porous silica	INERTCLONE ODS(3)	n/a	Reprosil Pur Basic C18	
PHENOMENEX	Fully-porous silica	JUPITER C4	n/a	ReproSil-XR 300 C4	Prosphere HP 300 C4
PHENOMENEX	Fully-porous silica	JUPITER C18	n/a	ReproSil-XR 300 C18	Prosphere HP 300 C18
PHENOMENEX	Fully-porous silica	JUPITER C8	n/a	ReproSil-XR 300 C8	n/a
PHENOMENEX	Fully-porous silica	JUPITER Proteo	n/a	Reprosher 100 C12	n/a
PHENOMENEX	Core-Shell	KINETEX Biphenyl	n/a	Reproshell Biphenyl	n/a
PHENOMENEX	Core-Shell	KINETEX C18	n/a	Reproshell ODS-3	Reproshell ODS-1
PHENOMENEX	Core-Shell	KINETEX C8	n/a	Reproshell C8	n/a
PHENOMENEX	Fully-porous silica	IB-Sil C18	n/a	Reprosil 120 XR C18	n/a
PHENOMENEX	Fully-porous silica	IB-Sil C18 BDS	n/a	ReproSil 120 BDS C18	n/a
PHENOMENEX	Fully-porous silica	IB-Sil C8	n/a	Reprosil 120 XR C8	n/a
PHENOMENEX	Fully-porous silica	IB-Sil C8 BDS	n/a	ReproSil 120 BDS C8	n/a
PHENOMENEX	Core-Shell	KINETEX HILIC	n/a	Reproshell HILIC	n/a
PHENOMENEX	Core-Shell	KINETEX PAH	n/a	n/a	
PHENOMENEX	Core-Shell	KINETEX PFP/F5	n/a	Reproshell PFP	n/a
PHENOMENEX	Core-Shell	KINETEX Phenyl-Hexyl	n/a	Reproshell Phenylhexyl	n/a
PHENOMENEX	Core-Shell	KINETEX Polar C18	n/a	n/a	n/a
PHENOMENEX	Core-Shell	KINETEX PS C18	n/a	n/a	n/a
PHENOMENEX	Core-Shell	KINETEX XB-C18	n/a	ReproShell SB-C18	Reproshell ODS-3
PHENOMENEX	Core-Shell	KINETEX EVO	n/a	Reprosil Phoenix	Reproshell ODS-3
PHENOMENEX	Fully-porous silica	LUNA C18	n/a	Reprosil Pur Basic C18	n/a
PHENOMENEX	Fully-porous silica	LUNA C18(2)	Reprosil Saphir 100 C18	Reprospher 100 C18-DE	n/a
PHENOMENEX	Fully-porous silica	LUNA C18(2)-HST	Reprosil Saphir 100 C18	Reprospher 100 C18-DE	n/a
PHENOMENEX	Fully-porous silica	LUNA C5	Reprosil Pur Basic C5	n/a	n/a
PHENOMENEX	Fully-porous silica	LUNA C8	Reprosil Pur Basic C8	n/a	n/a
PHENOMENEX	Fully-porous silica	LUNA C8(2)	Reprosil Saphir 100 C8	n/a	n/a
PHENOMENEX	Fully-porous silica	LUNA CN	n/a	n/a	Reprospher 100 CN
PHENOMENEX	Fully-porous silica	LUNA HILIC	n/a	n/a	Reprospher 100 Diol
PHENOMENEX	Fully-porous silica	LUNA NH2	n/a	n/a	Reprospher 100 NH2
PHENOMENEX	Fully-porous silica	LUNA Omega C18	n/a	Reprosil Saphir 100 C18	Reprosil XR 120 C18
PHENOMENEX	Fully-porous silica	LUNA Omega Polar	n/a	Reprosil XR 120 C18-MS	n/a
PHENOMENEX	Fully-porous silica	LUNA Omega PS	n/a	n/a	Reprospher 120 SAX C8
PHENOMENEX	Fully-porous silica	LUNA Omega Sugar	n/a	n/a	n/a
PHENOMENEX	Fully-porous silica	LUNA PFP(2)	n/a	Reprospher 100 PFP	n/a
PHENOMENEX	Fully-porous silica	LUNA Phe-Hexyl	Reprospher 100 Phenyl-hexyl-e	n/a	n/a
PHENOMENEX	Fully-porous silica	LUNA SCX	Reprosil Saphir 100 SCX	n/a	n/a

PHENOMENEX	Fully-porous silica	LUNA Si(2)	Reprosil Saphir 100 Si	n/a	n/a
PHENOMENEX	CHIRAL - AMYLOSE	LUX i-Amylose-1	Reprosil Chiral-MIA	n/a	n/a
PHENOMENEX	CHIRAL- CELLULOSE	LUX i-Cellulose-5	Reprosil Chiral-MIC	n/a	n/a
PHENOMENEX	CHIRAL - AMYLOSE	LUX Amylose-1	Reprosil Chiral-AM	n/a	n/a
PHENOMENEX	CHIRAL - AMYLOSE	LUX Amylose-2	Reprosil Chiral-YM	n/a	n/a
PHENOMENEX	CHIRAL- CELLULOSE	LUX Cellulose-3	Reprosil Chiral-JM	n/a	n/a
PHENOMENEX	CHIRAL- CELLULOSE	LUX Cellulose-1	Reprosil Chiral-OM	n/a	n/a
PHENOMENEX	CHIRAL- CELLULOSE	LUX Cellulose-4	Reprosil Chiral-XM	n/a	n/a
PHENOMENEX	CHIRAL- CELLULOSE	LUX Cellulose-2	Reprosil Chiral-ZM	n/a	n/a
PHENOMENEX	CHIRAL- CELLULOSE	LUX AMP	n/a	n/a	n/a
PHENOMENEX	MONOLITHIC SILICA	ONYX	n/a	n/a	n/a
PHENOMENEX	GPC	PHENOGEL 50A	ReproGel GPC, 50 A°	n/a	n/a
PHENOMENEX	GPC	PHENOGEL 100A	ReproGel GPC, 100 A°	n/a	n/a
PHENOMENEX	GPC	PHENOGEL 500A	ReproGel GPC, 500 A°	n/a	n/a
PHENOMENEX	GPC	PHENOGEL 1.000 A	ReproGel GPC, 1.000 A°	n/a	n/a
PHENOMENEX	GPC	PHENOGEL 10.000A	ReproGel GPC, 10.000 A°	n/a	n/a
PHENOMENEX	GPC	PHENOGEL 100.000A	ReproGel GPC, 100.000	n/a	n/a
PHENOMENEX	GPC	PHENOGEL 1.000.000A	ReproGel GPC, 1.000.000 A°	n/a	n/a
PHENOMENEX	GPC	PHENOGEL Linear	ReproGel GPC, 5 µm, linear (Mixed-bed)	n/a	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE-Next C18	n/a	Reprosil 80 ODS-1	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE-Next Si	n/a	Reprosil 80 Si	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE-Next Phenyl	n/a	Reprosil 80 Phenyl	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE-Next C8	n/a	Reprosil 80 C8	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE NH2	n/a	Reprosil 80 NH2	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE C1	n/a	Reprosil 80 C1	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE C6	n/a	Reprosil 80 Hexyl (C6)	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE SCX	n/a	Reprosil 80 SCX	n/a
PHENOMENEX	Fully-porous silica	PHENOSPHERE SAX	n/a	Reprosil 80 SAX	n/a
PHENOMENEX	PSDVB	POLYMERX	Reprogel 100 RP-1	Repromer 100 RPS	n/a
PHENOMENEX	Polymer	POLYSEP-GFC-P 1000	n/a	Repromer OH-2500	n/a
PHENOMENEX	Polymer	POLYSEP-GFC-P 2000	n/a	n/a	n/a
PHENOMENEX	Polymer	POLYSEP-GFC-P 3000	n/a	Repromer OH-3000	n/a
PHENOMENEX	Polymer	POLYSEP-GFC-P 4000	n/a	Repromer OH-4000	n/a
PHENOMENEX	Polymer	POLYSEP-GFC-P 5000	n/a	Repromer OH-5000	n/a
PHENOMENEX	Polymer	POLYSEP-GFC-P 6000	n/a	n/a	n/a
PHENOMENEX	Polymer	POLYSEP-GFC-P Linear	n/a	n/a	n/a
PHENOMENEX	Fully-porous silica	PRODIGY ODS(2)	n/a	Reprosil Pur 120 C18	n/a
PHENOMENEX	Fully-porous silica	PRODIGY C8	n/a	Reprosil Pur 120 C8	n/a
PHENOMENEX	Fully-porous silica	PRODIGY ODS(3)	n/a	Reprosil Pur Basic C18	n/a
PHENOMENEX	Fully-porous silica	PRODIGY Phenyl (PH-3)	n/a	ReproSil-Pur 120 Phenyl	n/a
PHENOMENEX	SUGAR	REZEX RAM Carbohydrate	Repromer LEX Ag	n/a	n/a
PHENOMENEX	SUGAR	REZEX RCM Monosaccharide	Repromer Ca	Reprogel Ca	n/a
PHENOMENEX	SUGAR	REZEX RCU Sugar Alcohols	Repromer Ca	Reprogel Ca	n/a
PHENOMENEX	SUGAR	REZEX RFQ Fast Acid	Repromer H	Reprogel H	n/a
PHENOMENEX	SUGAR	REZEX RHM Monosaccharide	Repromer H	Reprogel H	n/a
PHENOMENEX	SUGAR	REZEX RNM Carbohydrate	Repromer Na	Reprogel Na	n/a
PHENOMENEX	SUGAR	REZEX RNO Oligosaccharide	Repromer Na	Reprogel Na	n/a
PHENOMENEX	SUGAR	REZEX ROA Organic Acid	Repromer H	Reprogel H	n/a

PHENOMENEX	SUGAR	REZEX RPM Monosaccharide	Repromer Pb	Reprogel Ph	n/a
PHENOMENEX	SUGAR	REZEX RSO Oligosaccharide	Repromer LEX Ag	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHERECLONE Silica	Reprosil 80 Si	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHERECLONE C6	Reprosil 80 C6	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHERECLONE C8	Reprosil 80 C8	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHERECLONE ODS(1)	Reprosil 80 ODS-1	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHERECLONE ODS(2)	Reprosil 80 ODS-2	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHERECLONE NH2	Reprosil 80 NH2	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHEREX C8	n/a	Reprosil 100 C8	n/a
PHENOMENEX	Fully-porous silica	SPHERECLONE SAX	Reprosil 80 SAX	n/a	n/a
PHENOMENEX	Fully-porous silica	SYNERGI FUSION-RP	n/a	n/a	Platinum EPS C18
PHENOMENEX	Fully-porous silica	SYNERGI HYDRO-RP	n/a	Reprosil Pur 120 C18-aq	Reprospher 100 C18-aqua
PHENOMENEX	Fully-porous silica	SYNERGI MAX-RP	n/a	Reprospher 100 C12	n/a
PHENOMENEX	Fully-porous silica	SYNERGI POLAR-RP	n/a	Reprosil Star Phenyl-ether	n/a
PHENOMENEX	Fully-porous silica	ULTRACARB C8	Reprosil 80 C8	n/a	n/a
PHENOMENEX	Fully-porous silica	ULTRACARB ODS(20)	Equisorb ODS-20	n/a	n/a
PHENOMENEX	Fully-porous silica	ULTRACARB ODS(30)	Equisorb ODS-30	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHEREX C6	Reprosil 80 C6	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHEREX NH2	Reprosil 80 NH2	n/a	n/a
PHENOMENEX	Fully-porous silica	SPHEREX C18	Reprosil 80 C18	n/a	n/a
PHENOMENEX	Fully-porous silica	ULTREMEX C8	Reprosil 80 C8	n/a	n/a
PHENOMENEX	Fully-porous silica	ULTREMEX C18	Reprosil 80 C18	n/a	n/a
PHENOMENEX	SEC-Silica	YARRA SEC-200	Reprosil 200 SEC	n/a	n/a
PHENOMENEX	SEC-Silica	YARRA SEC-3000	Reprosil 300 SEC	n/a	n/a
THERMO	Fully-porous silica	Acclaim	n/a	n/a	Reprosil Pur Basic C18-HD
THERMO	Fully-porous silica	Acclaim surfactant	n/a	n/a	Ultrasep ES Tensid
ACE	Fully-porous silica	ACE C18 amid	n/a	Stability 100 Amid-C18	n/a
ACE	Fully-porous silica	ACE C18 AR	n/a	Reprospher 100 C18-Phenyl	n/a
ACE	Fully-porous silica	ACE C18 PFP	n/a	Exsil Plus C18 PFP	n/a
WATERS	Ethylene-bridged silica	Acquity BEH	n/a	Reprosil PHOENIX	n/a
AGILENT	Fully-porous silica	Agilent SB-AQ	n/a	n/a	Platinum C18
BIORAD	SUGAR	Aminex HPX-87C Column	Repromer Ca	n/a	n/a
BIORAD	SUGAR	Aminex HPX-87H Column	Repromer H	n/a	n/a
BIORAD	SUGAR	Aminex HPX-87K Column	Repromer K	n/a	n/a
BIORAD	SUGAR	Aminex HPX-87N Column	Repromer Na	n/a	n/a
BIORAD	SUGAR	Aminex HPX-87P Column	Repromer Pb	n/a	n/a
WATERS	Fully-porous silica	Atlantis T3	n/a	Reprosil Pur 120 C18-aq	Reprospher 100 C18-aqua
Waters	Fully-porous silica	Bondapak	Bondapak	Reprobond	n/a
MERCK /SUPELCO	Fully-porous silica	Cyclobond I 2000	n/a	n/a	Reprosil Chiral-Beta-CD
GRACE	Fully-porous silica	Econospher	Exsil 80	n/a	n/a
HAMILTON	Polymer	Hamilton PRP-1	n/a	Repromer 100 RPS	n/a
HAMILTON	Polymer	HAMILTON PRP-X200	n/a	n/a	Repromer SCX, 10 µm
THERMO	Fully-porous silica	Hypersil BDS	Hypersil BDS	Equisil BDS	Reprosil 120 BDS
THERMO	Fully-porous silica	Hypersil C18	Hypersil C18	Equisil ODS	n/a
THERMO	Fully-porous silica	Hypersil Gold	n/a	Equisil Gold	Equisil Gold
THERMO	Fully-porous silica	Hypersil green PAH	n/a	Ultrasep ES PAH Reprosil PAH	n/a
GL Science	Fully-porous silica	Inertsil ODS	Inertsil ODS	n/a	Reprosil-Pur 60 ODS-3

GL Science	Fully-porous silica	Inertsil ODS-2	Inertsil ODS-2	ReproSpher 100 C18-DE	ReproSil-Pur 120 ODS-3 ReproSil Saphir 100 C18
GL Science	Fully-porous silica	Inertsil ODS-3	Inertsil ODS-3	ReproSil-Pur Basic C18	ReproSil-Pur Basic-C18, 5 µm
NOURYON	Fully-porous silica	Kromasil	Kromasil	Reprokrom	n/a
MERCK /SUPELCO	Fully-porous silica	Lichrospher RP18	n/a	ReproSpher 100 C18-TN, 5 µm	n/a
MERCK /SUPELCO	Fully-porous silica	Lichrospher®100 RP-18 endcapped	n/a	ReproSpher 100 C18-TDE, 5 µm	n/a
Waters	Fully-porous silica	Nova-Pack	Repropak	n/a	n/a
MACHEREY-NAGEL	Fully-porous silica	Nucleosil 300	n/a	n/a	ReproSil-XR 300 C18
MACHEREY-NAGEL	Fully-porous silica	Nucleosil C18	n/a	ReproSpher 100 C18	n/a
WHATMAN	Fully-porous silica	Partisil C18	Partisil C18	Repropart ODS-3	n/a
WHATMAN	Fully-porous silica	Partisil SCX	Partisil SCX	Repropart SCX	n/a
WHATMAN	Fully-porous silica	Partisphere SCX	n/a	ReproSil XR 120 SCX	n/a
RESTEK	Fully-porous silica	Pinnacle (II) PAH	n/a	Ultrasep ES PAH Reprosil PAH	n/a
MERCK /SUPELCO	Fully-porous silica	Purospher Star RP18e	ReproSil Star RP18-e	n/a	n/a
Waters	Fully-porous silica	Spherisorb	Spherisorb	Reprosil 80	n/a
WATERS	Fully-porous silica	Sunfire	n/a	ReproSpher 100 DE	n/a
MERCK /SUPELCO	Fully-porous silica	Supelcosil ABZ, 200A°	n/a	n/a	ReproSil Amid-C18 ABZ
MERCK /SUPELCO	Fully-porous silica	SUPELCOSIL™ LC-8	n/a	Reprosil XR C8	Equisil
WATERS	Fully-porous silica	Symmetry	n/a	ReproSil-Pur Basic	n/a
WATERS	Fully-porous silica	Symmetry Shield	n/a	n/a	Stability 120 BS-C17
THERMO	Fully-porous silica	Thermo Scientific Acclaim 120 C18	n/a	n/a	Reprosil Pur 120 ODS-3
GRACE	Fully-porous silica	Vydac 201TP C18	n/a	n/a	ReproSpher 300 C18-TN
HICHROM	Fully-porous silica	Vydac 201TP C18	n/a	n/a	ReproSpher 300 C18-TN
GRACE	Fully-porous silica	Vydac 208TP D5C8	n/a	Prosphere 300 C8	Reprosil XR 300 C8
HICHROM	Fully-porous silica	Vydac 208TP D5C8	n/a	Prosphere 300 C8	Reprosil XR 300 C8
GRACE	Fully-porous silica	Vydac 214TP D5C4	n/a	Prosphere 300 C4	Reprosil XR 300 C4
HICHROM	Fully-porous silica	Vydac 214TP D5C4	n/a	Prosphere 300 C4	Reprosil XR 300 C4
GRACE	Fully-porous silica	Vydac 218TP C18	n/a	Prosphere 300 C18	Reprosil XR 300 C18
HICHROM	Fully-porous silica	Vydac 218TP C18	n/a	Prosphere 300 C18	Reprosil XR 300 C18
WATERS	Ethylene-bridged silica	Xbridge	n/a	Reprosil PHOENIX	n/a
WATERS	Fully-porous silica	XSelect C18 CSH	n/a	n/a	Avanti BDS Quattro
WATERS	Ethylene-bridged silica	Xterra	n/a	Reprosil PHOENIX	n/a
YMC	Fully-porous silica	YMC Pack ODS-A	ReproSil-Pur 120 ODS-3	n/a	n/a
YMC	Fully-porous silica	YMC-Pack ProC18	n/a	ReproSil 100 C18-XBD	ReproSil 100 C18-XBD
AGILENT	Fully-porous silica	Zorbax 300 SCX	n/a	ReproSil XR 300 SCX	Exsil Plus 300 SCX
AGILENT	Fully-porous silica	Zorbax Bonus RP	n/a	n/a	Stability 120 Amid-C12
AGILENT	Fully-porous silica	Zorbax Eclipse PAH	n/a	Ultrasep ES PAH Reprosil PAH	n/a
AGILENT	Fully-porous silica	Zorbax SB-C18	n/a	Reprosil Star SB-C18	n/a
THERMO	Core-Shell	Accucore	n/a	Reproshell	n/a
GL Science	Fully-porous silica	Inertsil ODS-HL – C18	n/a	n/a	ReproSil-Pur Basic-C18-HD
GL Science	Fully-porous silica	Inertsil ODS-SP – C18	n/a	n/a	Reprosil Pur 120 C8